


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Etidronic acid - 60% aqueous solution

Product Number: **H6773**
CAS Number: 2809-21-4
MDL: MFCD00070585
Formula: C₂H₈O₇P₂
Formula Weight: 206.03 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Colorless to Faint Yellow
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Phosphorus (P)	16 - 20 %
Purity (Titration by NaOH)	55.0 - 65.0 %
Recommended Retest Period	-----
6 Years	

Specification: PRD.0.ZQ5.10000025346

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.