

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
cis-3-Hexen-1-ol - ≥98%, FCC, FG

Product Number: W256307
CAS Number: 928-96-1
MDL: MFCD00063217
Formula: C6H12O
Formula Weight: 100.16 g/mol



TEST	Specification
Appearance (Color)	Colorless to Faint Yellow
Appearance (Form)	Liquid
Refractive index at 20 °C	1.436 - 1.443
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 98.0 %
Sum Of Isomers Cis and Trans	
GC (area %)	≥ 92.0 %
Cis Isomer	
Specific Gravity	0.846 - 0.850
At 25°C	
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period	-----
5 Years	

Specification: PRD.1.ZQ5.10000048833

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.