

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
NFκB (p65) Control Peptide human - ~1 mg/mL, buffered aqueous solution

Product Number: N6659
MDL: MFCD03457712
Storage Temperature: -20 °C

TEST	Specification
Formulation	-----
Solution containing 0.01 % Sodium Azide	
Protein Concentration	Conforms
Approx 1 mg/mL	
Application	Conforms
to block specific interaction of anti-NFκB p65 with the NFκB p65 subunit; should be used at 1 ug per 1 uL of antiserum used	
Specification: PRD.0.ZQ5.10000046758	

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.