

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Gallium – 99.9995% trace metals basis

**Product Number:****203319**

CAS Number:

7440-55-3

MDL:

MFCD00134045

Formula:

Ga

Formula Weight:

69.72 g/mol

Storage Temperature:

2 - 8 °C

Ga

TEST	Specification
------	---------------

Appearance (Color)

Conforms to Requirements

Silver or Grey

Appearance (Form)

Conforms to Requirements

Metallic Solid and/or Liquid

ICP Major Analysis

Confirmed

Confirms Gallium Component

Purity

Meets Requirements

99.9995% Based On Trace Metals Analysis

Trace Metal Analysis

≤ 5.5 ppm

Specification: PRD.0.ZQ5.10000077267

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.