

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti-Inositol 1,4,5-Triphosphate Receptor (Type I) antibody produced in rabbit – whole antiserum, buffered aqueous solution

Product Number:

I157

MDL:

MFCD01092546

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Formulation

Conforms

Phosphate buffered saline containing 0.05% sodium azide

Working Dilution

Conforms

By Direct Immunoblotting; assay dependent

Specification: PRD.1.ZQ5.10000054764

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.