

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tris Buffered Saline - 10 x, solution

Product Number: T5912

MDL: MFCD00132476

Formula Weight: 18.02 g/mol

Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Liquid
pH	7.2 - 7.6
Conductivity	100 - 130 mS/cm
Dot Blot	≤ 2.5 ng

Specification: PRD.7.ZQ5.10000006660

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.