

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

PSS-(2-(trans-3,4-Cyclohexanediol)ethyl)-Heptaisobutyl substituted

Product Number:

560286

CAS Number:

480439-48-3

MDL:


MFCD03427218

Formula:

C₃₆H₇₈O₁₄Si₈

Formula Weight:

959.68 g/mol



TEST

Specification

Appearance (Color)

White to Off-White

Appearance (Form)

Powder

Melting Point

156.0 - 165.0 °C

Infrared spectrum

Conforms to Structure

Proton NMR spectra

Conforms to Structure

Carbon

42.3 - 47.7 %

Specification: PRD.0.ZQ5.10000028274

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.