


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Mefenamic acid

Product Number: **M4267**
CAS Number: 61-68-7
MDL: MFCD00051721
Formula: C₁₅H₁₅NO₂
Formula Weight: 241.29 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Titration with NaOH	≥ 98
Purity (HPLC)	≥ 98 %
Solubility (Color)	Faint Yellow
Solubility (Turbidity)	Clear to Slightly Hazy
50 mg/ml, 1 M Ammonium Hydroxide	

Specification: PRD.1.ZQ5.10000045649

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.