

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

LB Agar Ampicillin - 50, Plates - pre-poured agar plates with 50 µg/ml ampicillin

Product Number: **L0168**

Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance Amber agar plates	Conforms
Performance Testing Growth promotion and Antibiotic or Amino acid sensitivity	Pass
Sterility	Pass

Specification: PRD.0.ZQ5.10000037327

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.