

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Phosphatase Inhibitor Cocktail 3 - DMSO solution

Product Number: **P0044**

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Form) Clear Solution	Conforms
Appearance (Colour)	Colorless to Faint Yellow
Enzym. Assay (Name)	Alkal. Phosphatase or PP1-Alpha
Enzym. Assay (Inhib. Activity) (With 1% inhibitor)	≥ 50.00 %
Recommended Retest Period 2 Years	-----

Specification: PRE.2.ZQ5.10000001024

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.