

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Estriol-2,3,4-¹³C₃ - 99 atom % ¹³C, 97% (CP)

Product Number: **731668**
CAS Number: 1255639-56-5

Formula: 13C₃C₁₅H₂₄O₃
Formula Weight: 291.33 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Purity (HPLC)	≥ 97.0 %
Proton NMR Spectrum	Conforms
Carbon NMR Spectrum	Conforms
LCMS	Conforms
¹³ C Enrichment	≥ 99.0 %

Specification: PRD.0.ZQ5.10000102097

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.