


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Bis[2-(methacryloyloxy)ethyl] phosphate

Product Number: 496758
CAS Number: 32435-46-4
MDL: MFCD00128922
Formula: C₁₂H₁₉O₈P
Formula Weight: 322.25 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)
may contain black particles
Appearance (Form)
Infrared spectrum
Proton NMR spectrum

Colorless to Yellow
Liquid
Conforms to Structure
Conforms to Structure

Specification: PRD.1.ZQ5.1000005424

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.