

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Activated charcoal – DARCO®, 20–40 mesh particle size, granular

**Product Number:** 242268  
**CAS Number:** 7440-44-0  
**MDL:** MFCD00133992  
**Formula:** C  
**Formula Weight:** 12.01 g/mol

C

TEST	Specification
Appearance (Color)	Black
Appearance (Form)	Crystal or Granules
Moisture Content	≤ 12 %
Dust	≤ 0.4 %
Particle Size Distribution	≤ 10 %
+ 20 Mesh	
Particle Size Distribution	≤ 5 %
-40 Mesh	
Note	Conforms
Darco is a Registered Trademark of Amercian Norit Company, Inc.	

Specification: PRD.0.ZQ5.10000024857

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.