

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

N,N,N',N'-Tetramethylethylenediamine – BioReagent, suitable for electrophoresis, ~99%

Product Number:

T9281

CAS Number:

110-18-9

MDL:


MFCD00008335

Formula:

C₆H₁₆N₂

Formula Weight:

116.2 g/mol



TEST

Specification

Appearance (Clarity)

Clear

Appearance (Color)

Colorless to Faint Yellow

Appearance (Form)

Liquid

Suitability

Suitable

Suitable for Use in Polyacrylamide Gel Electrophoresis

Purity (GC)

≥ 99 %

Specification: PRD.1.ZQ5.10000003730

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.