

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Phosphatase, Alkaline from calf intestine - buffered aqueous glycerol solution

Product Number: **P4978**

MDL: MFCD00131849

Storage Temperature: -20 °C

TEST	Specification
units/ml	≥ 10000
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Liquid
DNase, Exonuclease Detection Using 5 units/reaction when incubated at 37 degrees Celcius for 1 hour.	None Detected
RNAse Detection Using 5 units/reaction when incubated at 37 degrees Celcius for 1 hour.	None Detected
NICKase, Endonuclease Detection Using 5 units/reaction when incubated at 37 degrees Celcius for 1 hour.	None Detected

Specification: PRD.3.ZQ5.10000005514

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.