

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com

Email USA: techserv@sial.com


Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Spiro-MeOTAD - 99% (HPLC)

Product Number: **792071**
CAS Number: 207739-72-8

Formula: C₈₁H₆₈N₄O₈
Formula Weight: 1,225.43 g/mol



TEST	Specification
Appearance (Color)	White to Faint Green to Faint Yellow
Appearance (Form)	Powder or Crystals
Proton NMR spectra	Conforms to Structure
Purity (HPLC)	≥ 98.5 %
Carbon	75.0 - 83.8 %
Nitrogen	
Specification: Report Results	
Wavelength	290 - 310 nm
UV/Vis Absorbance Max	
0.01 g/L CH ₂ Cl ₂	
Wavelength	370 - 390 nm
UV/Vis Absorbance Max	
Recommended Retest Period	-----
2 years	

Specification: PRD.3.ZQ5.10000047904

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.