


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: H3128
CAS Number: 3737-95-9
MDL: MFCD00004078
Formula: C₂₁H₁₄N₂O₇S
Formula Weight: 438.41 g/mol



TEST	Specification
Appearance (Color)	Black
Appearance (Form)	Powder
Solubility (Color) 10 mg/ml, 1 M NaOH	Blue to Dark Blue
Suitability Suitable for use as an indicator for the titration of calcium.	Suitable

Specification: PRD.0.ZQ5.10000050605

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.