

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Chemiluminescent Peroxidase Substrate for ELISA – sufficient for 60 mL substrate

Product Number:

CPS260

Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

Kit Comp (Component 1)

Conforms

Apperance: Clear colorless solution

pH: 7.25-7.75

Kit Comp (Component 2)

Conforms

Appearance: Clear colorless to pinkish tint liquid

pH: 9.35-9.85

Specification: PRD.4.ZQ5.10000006470

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.