


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Heptane - ReagentPlus[®], 99%

Product Number: H2198
CAS Number: 142-82-5
MDL: MFCD00009544
Formula: C₇H₁₆
Formula Weight: 100.2 g/mol



| TEST | Specification |
|---------------------------|-----------------------|
| Appearance (Color) | Colorless |
| Appearance (Form) | Liquid |
| Infrared spectrum | Conforms to Structure |
| Refractive index at 20 °C | 1.385 - 1.389 |
| Purity (GC) | ≥ 98.5 % |

Specification: PRD.1.ZQ5.10000003992

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.