

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
DBE-3

Product Number: 422037

MDL: MFCD00191969



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
¹ H NMR Spectrum	Conforms to Structure
APHA Color	≤ 15 APHA
Dimethyl Glutarate	5 - 15 %
Dimethyl Adipate	85 - 95 %
Dimethyl Succinate	≤ 1.0 %
Ester Content	≥ 99.0 %
Total Esters	
Product of Supplier	Confirmed
INVISTA	

Specification: PRD.0.ZQ5.10000099592

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.