

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

L-Glutathione reduced - ≥98.0%

Product Number: G4251 CAS Number: 70-18-8

 MDL:
 MFCD00065939

 Formula:
 C10H17N3O6S

 Formula Weight:
 307.32 g/mol

 Storage Temperature:
 2 - 8 ℃

TEST

Appearance (Color) Appearance (Form) Solubility (Color) Solubility (Turbidity) 50 mg/mL, H2O 1H NMR Spectrum

Purity (HPLC)
Titration with lodine

Recommended Retest Period

5 years

Specification: PRD.2.ZQ5.10000014977

Specification

White Powder Colorless Clear

Conforms to Structure

> 98.0 % > 99 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.