

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4,4'-Methylenebis(phenyl isocyanate) - 98%

Product Number: 256439
CAS Number: 101-68-8
MDL: MFCD00036131
Formula: C₁₅H₁₀N₂O₂
Formula Weight: 250.25 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form) Flakes or Chips	Conforms to Requirements
Proton NMR Spectrum	Conforms to Structure
Purity (HPLC)	≥ 97.5 %
Carbon	70.2 - 73.8 %
Nitrogen	10.9 - 11.5 %

Specification: PRD.4.ZQ5.10000003616

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.