

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Opoponax oil - naturalProduct Number: **W522600**
CAS Number: 8021-36-1
MDL: MFCD00678435

TEST	Specification
------	---------------

Appearance (Color)	Yellow to Very Dark Yellow to Very Dark Orange
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.485 - 1.595
Specific Gravity at 20 deg C	0.980 - 1.015
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period 5 Years	-----

Specification: PRD.4.ZQ5.10000041780

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.