


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
γ-Hexalactone - natural, 97%, FGProduct Number: **W255611**

CAS Number: 695-06-7

Formula: C₆H₁₀O₂
Formula Weight: 114.14 g/mol

TEST	Specification
------	---------------

Appearance (Color)	Colorless to Faint Yellow
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.434 - 1.442
Infrared spectrum	Complete
Purity (GC)	≥ 96.5 %
Optical Rotation	
C= Neat	
Specification: Report Result	
Arsenic (As)	≤ 3 ppm
Mercury (Hg)	≤ 1 ppm
Cadmium (Cd)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period	-----
60 Months	

Specification: PRD.4.ZQ5.10000034278

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.