

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Fumonisin B₂ from *Fusarium moniliforme* – sphingosine N-acyltransferase inhibitor

Product Number:

F3771

CAS Number:

116355-84-1

Formula:


C₃₄H₅₉NO₁₄

Formula Weight:

705.83 g/mol

Storage Temperature:

-20 °C

**TEST****Specification**

Appearance (Color)

White to Yellow to Orange

Appearance (Form)

Powder

Purity (TLC)

≥ 96 %

Solubility (Color)

Colorless to Faint Yellow

Solubility (Turbidity)

Clear

10 mg/ml, MeOH

Specification: PRD.0.ZQ5.10000101235

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.