

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Stemline® Keratinocyte Medium II - Serum-free, Without L-glutamine, calcium-free, with HEPES, Sodium Bicarbonate, and phenol red, liquid, sterile-filtered

Product Number: **S0196**

MDL: MFCD02098338

Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance Clear Pink/Brown	Conforms
pH	7.3 - 7.5
Sterility (0.2 micron filtration) sterile	Conforms
Growth promotion test Confirms	Conforms
Expiration Dated Product	-----

Specification: PRD.3.ZQ5.10000007534

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.