

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Dihydroceramide C2 - ≥97% (TLC), solid


Product Number: C7980

MDL: MFCD00236495

Formula: C₂₀H₄₁NO₃

Formula Weight: 343.54 g/mol

Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
5 mg/mL, EtOH	
Proton NMR spectrum	Conforms to Structure
Purity (TLC)	≥ 97 %

Specification: PRD.1.ZQ5.10000047912

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.