

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Monoclonal Anti-Insulin-like Growth Factor-I antibody produced in hamster - clone 126002, purified immunoglobulin, lyophilized powder

Product Number: **I9409**

Storage Temperature: -20 °C

TEST	Specification
Appearance (Form) from a 0.2 um filtered solution in PBS with 5% trehalose	Lyophilized Powder
Endotoxin Assay	Less than 0.1 Eu per 1 ug
Western Blot	1 to 2 ug/ml
Elisa Capture	2 to 8 ug/ml
Neutralization of Bioactivity	0.3 to 1.5 ug/ml

Specification: PRD.0.ZQ5.10000051548

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.