

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com


Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Folinic acid calcium salt hydrate – analytical standard

Product Number: **F7878**

MDL: MFCD00006704
Formula: C₂₀H₂₁CaN₇O₇ · xH₂O
Formula Weight: 511.50 g/mol



TEST	Specification
Appearance (Color)	Yellow to Orange
Appearance (Form)	Powder
Water (by Karl Fischer)	≤ 17.0 %
Solubility (Color)	Yellow to Green
Solubility (Turbidity)	Clear to Slightly Hazy
100 mg/mL, H ₂ O	
Calcium (Ca) anhydrous	7 - 9 %
Purity (HPLC)	≥ 90 %
Recommended Retest Period	-----
3 years	

Specification: PRD.0.ZQ5.10000069329

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.