

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Protein A-Sepharose® 4B, Fast Flow from Staphylococcus aureus – aqueous ethanol suspension

Product Number: P9424
MDL: MFCD00165563
Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color)	White
Appearance (Form)	Suspension
ml suspension/ml packed gel	1.0 - 1.4
mg/ml binding capacity	≥ 30
Human IgG	
Registered Trademark	Confirmed
Sepharose is a Registered Trademark of GE Healthcare	
Recommended Retest Period	-----
5 years	
Specification: PRD.0.ZQ5.10000027179	

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.