

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Poly(methyl methacrylate) – analytical standard, for GPC, average M_w 97,000 (Typical), average M_n 46,000 (Typical)

Product Number:

370037

CAS Number:


9011-14-7

MDL:

MFCD00134349

Formula:

$[\text{CH}_2\text{C}(\text{CH}_3)(\text{CO}_2\text{CH}_3)]_n$



TEST	Specification
Appearance (Color)	White
Appearance (Form) Crystals or Crystalline Powder	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Molecular Weight Typically 97,000	
Molecular Number Typically 46,000	
Note Lot Specific Analysis Included With Each Unit	Confirmed

Specification: PRD.2.ZQ5.10000033853

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.