

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

Benzophenone purified by sublimation, $\geq 99\%$

Product Number:

427551

CAS Number:

119-61-9

MDL:


MFCD00003076

Formula:

C₁₃H₁₀O

Formula Weight:

182.22 g/mol

**TEST****Specification**

Appearance (Color)

White

Appearance (Form)

Conforms to Requirements

Crystals or Powder

Infrared spectrum

Conforms to Structure

Purity (GC)

 $\geq 99.0\%$ **Remarks:**

Specification Date : 12/15/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.