


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Bromo-3-methylpyridine - 95%

Product Number: 448109
CAS Number: 3430-17-9
MDL: MFCD00239380
Formula: C₆H₆BrN
Formula Weight: 172.02 g/mol



TEST

Specification

Appearance (Color)	Colorless to Yellow to Brown
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Proton NMR spectra	Conforms to Structure
Titration by AgNO ₃	94.0 - 106.0 %
Purity (GC)	≥ 94.0 %

Specification: PRD.0.ZQ5.10000057127

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.