

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Nde I from *Neisseria denitrificans* – Restriction Enzyme

Product Number:

R5509

CAS Number:

84628-87-5

MDL:

MFCD00166615

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Enzymatic Activity (Un/mL)

≥ 10000

Endonuclease

Conforms

Up to 20U with lambda-DNA Not detectable

Up to 20U with pBR 322-DNA Not detectable

Exonuclease

Conforms

Up to 30U Not detectable

Specification: PRD.1.ZQ5.10000030615

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.