

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: 905178



Formula: C₃₉H₆₁N₇O₆S
Formula Weight: 756.02 g/mol
Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color)	White to Yellow
Appearance (Form)	Solid or Semi-solid
¹ H NMR Spectrum	Conforms to Structure
Carbon Content	58.9 - 65.1 %
Nitrogen Content	12.3 - 13.6 %
Sulfur Content	
Specification: Report Only	
Specific Rotation	
0.1%; DMF	
Specification: Report Only	
LCMS	Confirmed
Confirms M+	
Purity (HPLC)	≥ 95.0 %

Specification: PRD.0.ZQ5.10000101388

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.