

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

1,1,1,5,5,6,6,6-Octafluoro-2,4-hexanedione - technical grade, 90%

Product Number:

429503

CAS Number:

20825-07-4

MDL:


MFCD00041014

Formula:

C₆H₂F₈O₂

Formula Weight:

258.07 g/mol



TEST

Specification

Appearance (Color)

Colorless to Light Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Assay

≥ 89.5 %

Specification: PRD.2.ZQ5.10000078209

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

