


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium sulfate - ≥99.99% trace metals basis

Product Number: 204447
CAS Number: 7757-82-6
MDL: MFCD00003504
Formula: Na₂O₄S
Formula Weight: 142.04 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Chunks
Gravimetric Analysis % SO ₄ with Barium Chloride	66.3 - 69.0 %
ICP Major Analysis Confirms Na Component	Confirmed
Trace Metal Analysis	≤ 100.0 ppm
Purity ≥99.99% Based on Trace Metals Analysis	Meets Requirements

Specification: PRD.0.ZQ5.10000027899

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.