


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,3-Bis(tributylstannyl)benzene

Product Number: 721190
CAS Number: 24344-61-4
Formula: C₃₀H₅₈Sn₂
Formula Weight: 656.20 g/mol



TEST	Specification
Appearance (Color)	Colorless to Yellow
Appearance (Form)	Liquid
Proton NMR spectrum	Conforms to Structure
ICP Major Analysis % Sn	31.9 - 40.5 %

Specification: PRD.0.ZQ5.10000054658

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.