

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Bmf, N-Terminal antibody produced in rabbit – affinity isolated antibody, buffered aqueous solution

Product Number: **B1684**

MDL: MFCD04040203

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Formulation	Conforms
1 mg/ml Solution in phosphate buffered saline with 0.02% sodium azide	
Working Dilution	Conforms
by Direct Immunoblotting	
2 ug/mL using human HEPG2 Hepatoblastoma and 293 Embryonic kidney cell lysates	

Specification: PRD.0.ZQ5.10000091204

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.