

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Silver hexafluorophosphate - 99.99% trace metals basis

Product Number: 450219
CAS Number: 26042-63-7
MDL: MFCD00003415
Formula: AgF₆P
Formula Weight: 252.83 g/mol



TEST	Specification
Appearance (Color)	Off-White
Appearance (Form)	Crystalline Powder
Titration by KSCN % Ag	41.8 - 43.6 %
ICP Major Analysis Confirms Silver and Phosphorus Component	Confirmed
Purity 99.99% Based On Trace Metals Analysis	Meets Requirements
Trace Metal Analysis	≤ 150.0 ppm

Specification: PRD.0.ZQ5.10000066094

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.