

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti-Mouse IgG (whole molecule)-FITC antibody produced in goat - IgG fraction of antiserum, buffered aqueous solution

Product Number: **F9006**

MDL: MFCD00162642

Storage Temperature: -20 °C

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Faint Yellow to Yellow
Appearance (Form)	Solution
F/P Molar Ratio	3.0 - 5.0
A280/A496 UV Absorption	1.0 - 1.5
Protein Content	10 - 20 mg/ml
Immunofluorescence (Direct)	Minimum Dilution 1:64
Using Mouse Spleen Cells	

Specification: PRD.3.ZQ5.10000007733

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.