

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


Potassium 6-chlorohexanoyltrifluoroborate - ≥95%

Product Number:

900609

CAS Number:

1370291-67-0



Formula:

C₆H₁₀BClF₃KO

Formula Weight:

240.50 g/mol

Storage Temperature:

2 - 8 °C

TEST

Specification

Appearance (Color)

White to Off White

Appearance (Form)

Powder

1H NMR Spectrum

Conforms to Structure

11B NMR Spectrum

Conforms to Structure

19F NMR Spectrum

Conforms to Structure

Carbon Content

28.5 - 31.5 %

Specification: PRD.0.ZQ5.10000083157

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.