


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Octanoic acid - ≥99%

Product Number: C2875
CAS Number: 124-07-2
MDL: MFCD00004429
Formula: C₈H₁₆O₂
Formula Weight: 144.21 g/mol



TEST	Specification
Appearance (Color)	Colorless to Faint Yellow
Appearance (Turbidity)	Clear
Appearance (Form)	Liquid
Proton NMR Spectrum	Conforms to Structure
GC (area %)	≥ 99.0 %
Recommended Retest Period	-----
4 Years	

Specification: PRD.1.ZQ5.10000019065

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.