

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Lithium chloride - anhydrous, beads, -10 mesh, 99.998% trace metals basis

Product Number: 429457
CAS Number: 7447-41-8
MDL: MFCD00011078
Formula: CLi
Formula Weight: 42.39 g/mol

LiCl

TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Beads
ICP Major Analysis	Confirmed
Confirms Lithium Component	
Purity	Conforms
99.998% Based On Trace Metals Analysis	
Trace Metal Analysis	≤ 25.0 ppm

Specification: PRD.0.ZQ5.10000084437

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.