

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
JumpStart™ REDTaq® ReadyMix™ Reaction Mix – for PCR

Product Number: P0982

Storage Temperature: -20 °C

TEST	Specification
Magnesium (Mg) by ICP	3.50 - 4.50 mM
DNAse, Exonuclease Detection	None Detected
RNAse Detection	None Detected
Suitable for PCR	Suitable
JumpStart Taq Antibody	Pass
Prevents non-specific product formation	

Specification: PRD.2.ZQ5.10000005967

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.