

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Germanium(II) bromide - 97%

Product Number: 572659
CAS Number: 24415-00-7
MDL: MFCD00049459
Formula: Br₂Ge
Formula Weight: 232.45 g/mol



TEST	Specification
Appearance (Color)	Yellow
Appearance (Form)	Powder
Titration with AgNO ₃ % Bromine	66.3 - 71.2 %
X-Ray Diffraction	Conforms to Structure

Specification: PRD.1.ZQ5.1000065061

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.