


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Phenol Red sodium salt

Product Number: P4758
CAS Number: 34487-61-1
MDL: MFCD00066901
Formula: C₁₉H₁₃NaO₅S
Formula Weight: 376.36 g/mol



TEST	Specification
Appearance (Color)	Red to Brown
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Solubility (Color)	Dark Red
Solubility (Form)	Solution
1 mg/mL, H ₂ O	
pH	Pass
6.8 (Yellow) to 8.2 (Red)	
Recommended Retest Period	-----
5 Years	

Specification: PRD.0.ZQ5.10000043318

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.