

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
nitro-Grela-I2 - ≥95%

Product Number: **901787**
CAS Number: 1874264-99-9

Formula: C₃₁H₃₇I₂N₃O₃Ru
Formula Weight: 854.52 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Brown to Very Dark Brown
Appearance (Form)	Powder
¹ H NMR Spectrum	Conforms to Structure
Titration with AgNO ₃ After Oxygen Combustion	95.0 - 105.0 %
ICP Major Analysis Confirms Ruthenium Component	Confirmed

Specification: PRD.0.ZQ5.10000102532

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

