


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Fluororesorcinol - 97%

Product Number: 679836
CAS Number: 103068-41-3
MDL: MFCD03789074
Formula: C₆H₅FO₂
Formula Weight: 128.10 g/mol
Storage Temperature: -20 °C



TEST	Specification
------	---------------

Appearance (Color)	White to Yellow and Faint Beige to Beige and Faint Brown to Brown and Faint Grey to Grey
Appearance (Form)	Solid or Powder or Crystals or Chunks
Infrared spectrum	Conforms to Structure
Carbon	54.3 - 58.2 %
GC (area %) (VPCT)	≥ 96.5 %

Specification: PRD.1.ZQ5.10000003623

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.