

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: **I7284**

Storage Temperature: -20 °C

TEST**Specification**

Appearance (Form)

from a 0.2 um filtered solution in PBS with 5% trehalose

Endotoxin Assay

Western Blot

Flow Cytometry

3-10 ug/ml per 1 million cells

Lyophilized Powder

Less than 0.1 Eu per 1 ug

0.1 to 0.2 ug/ml

Pass

Specification: PRD.0.ZQ5.10000100570

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.